

**SUBJECT:** Content Uniformity of 3.5mg Aquatabs®

**DATE:** 8<sup>th</sup> December, 2011


**FROM:** Ulick Stafford

**TO:** Whoever it may concern

**REFERENCE NUMBER:** MTS 010

Aquatabs® are manufactured to the highest quality standards in the industry. All tablets are made to GMP standards in our ISO 9000 facility. The evidence of our emphasis on quality can be measured in our products. We have recently developed a test of content uniformity of 3.5 mg effervescent Troclosen sodium tablets. The Content Uniformity of three random batches of Aquatabs® 3.5 mg were tested according to European Pharmacopeia 7.0. All tests complied with the standard and every sample was within 85% and 115 % of the average contents

For comparison we tested a sample of 3.5 mg effervescent Troclosen sodium tablets from another European manufacturer according to European Pharmacopeia 7.0. All samples were outside 75% - 125% specification for Content Uniformity.<sup>1</sup>



Signed by:

Ulick Stafford PhD, R&D Manager

on behalf of Medentech, Clonard Road, Wexford, Ireland.

<sup>1</sup> Medentech internal report, RD-CP-0024-Rev 1.